



Identification tool for vegetable and marine oils. (example is a marine oil).

ELUTION ORDER AND IDENTIFICATION ON
25M x 0,25mmID 0,2um film of polyethyleneglycol (wax column).


Errors may exist in the area from C16:1 to C16:4.




Complete chromatogram.




C12:0 to C16:1 section.




C16:0 – C18:1 section. (small peaks in the middle is a challenging area, may be deviations)



C18:0 – C20:1 section.



C20-C22:1 Section



C20:5 - C24:1 Section.